# FeLV Ag + FIV Ab combined test

Rev.: KI-40-10010-01

10 x #(2)800006

## FeLV Ag + FIV Ab combined test

Rapid test

- (PT) Teste FeLV Ag + FIV A
- (CS) Test FeLV Ag + FIV Ab
- (PL) Test FeLV Ag + FIV Ab
- (RO) FeLV Ag + FIV Ab Test
- (SK) Test FeLV Ag + FIV Ab
- (DA) FeLV Ag + FIV Ab Test
- (NO) FeLV Ag + FIV Ab Test
- (SV) FeLV Ag + FIV Ab-test
- (ZH) FeLV抗原+ FIV抗体测试

LOT

40100TTMMJJ

2022-03

30°C  
2°C

10x

10x

1x

1x

Distributed by:  
Covetrus BV  
Beversestraat 23  
5431 SL Cuijk (NL)  
cbproducts@covetrus.com

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
Rev.: PO-40-100-01

covetrus

## FeLV Ag + FIV Ab combined test double teststrip

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### (EN) FeLV Ag + FIV Ab combined test

#### Instructions for use

The Covetrus FeLV Ag + FIV Ab combined test is used to test for feline leucosis virus (FeLV) antigens and feline immunodeficiency virus (FIV) antibodies in serum, plasma or whole blood.

#### Handling of the sample

Do not keep the samples at room temperature for extended periods. Serum should be frozen at below -20°C for up to 3 days. For longer storage, the samples should be frozen at below -20°C. Whole blood should be stored between 2-8°C if the test is due to be carried out within two days from sampling. Bring the samples to room temperature before testing. Frozen samples must be fully thawed before testing.

General instructions for the use of blood samples:

• Separate the serum or plasma as rapidly as possible from the blood, to avoid haemolysis.

• Heparin blood and EDTA blood can be used to extract plasma. • Use only clear rather than haemolysed samples. • Testing should take place immediately after sampling.

#### Test procedure

##### Figure 1: Test de sérum et plasma

a. Draw a small sample from each tube on the serum or plasma. Take the test pipette and draw sample material up into pipette. Add one drop of the sample material to the two sample wells. Wait a few seconds; to allow the strip to fully absorb the sample material.

b. Then, open the reagent vessel and add two drops of reagent to each of the sample wells until the fluid begins to run up the test strip. If the liquid does not start to run up on the test strip after a few seconds, add another drop of reagent to the sample well.

**Figure 2: Draw of whole blood**

a. Draw blood from each tube on the test pipette and draw sample material up into pipette. Add one drop of the sample material to the two sample wells. Wait up to one minute, to allow the strip to fully absorb the sample material.

b. Then, open the reagent vessel and add two drops to each of the sample wells. The fluid should begin to run up the test strip. If the liquid does not start to run up on the test strip after a few seconds, add another drop of reagent to the sample well.

#### Analysis (Figures 3 and 4)

After 10 minutes, one or two red lines will appear in the reaction field, and the test can be analysed.

#### Positive test result (Figure 3 +):

Both test line and the control line appear in the reaction field.

The figure shows a positive FeLV antigen test result.

NB: Even a faint test line must be interpreted as a positive test result.

#### Negative test result (Figure 3 -):

Only the control line appears in the reaction field.

The figure shows a negative FeLV test result.

#### Invalid:

If no control line appears, test is invalid.

NB: All test results that are read after 10 minutes are invalid.

#### Disposal

All materials used must be disposed of carefully, because the material could be infectious.

Test materials and used cassettes must be put in plastic bags and then disposed of.

#### Plastic bags

For veterinary use only. • Only for personal use. • Use the cassette as soon as possible after opening. • Please follow the procedure for required sample quantity. • Do not apply any sample material to the reaction field. • Do not touch the reaction field. • Use only the fluid supplied. • Do not use after expiry date. • Do not use the test cassette if the packaging is damaged. • Observe analysis times indicated. • Store and transport the test only at the temperatures indicated.

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