Feline Leukemia Virus Antigen Test Kit

For veterinary use only.

English version

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Esclusivo per uso veterinario

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SNAP® FeLV

The SNAP FeLV Antigen test is an enzyme-linked immunosorbent assay (ELISA) for the detection of feline leukemia virus (FeLV) antigen in feline serum, plasma or whole blood. The detection of the p27 group-specific viral antigen (p27) is diagnostic for FeLV infection. The SNAP FeLV Antigen assay utilizes monoclonal antibodies to p27 and positive and negative controls.

The conjugate contains enzyme-conjugated antibody to p27. Upon mixing the conjugate and the test sample, conjugated monoclonal antibody will bind p27 antigen (if present). The sample/conjugate mixture is then added to the SNAP device and flows across the spotted matrix. The matrix-bound p27 antibody (FeLV Ag) will capture the p27-conjugated antibody complex. The device is then activated, rehydrating wash and substrate reagents stored within the device.

Color development in the FeLV Ag sample spot indicates the presence of FeLV antigen.

Precautions and warnings

- Use a separate sample tube and transfer pipet for each test.
- Using multiple transfer pipets and reagent bottles are not for quality control purposes. This may not interfere with the test results or interpretation.
- Samples, sample tubes, pipettes, conjugate and SNAP devices should be handled as though capable of transmitting FeLV. All waste should be properly decontaminated prior to disposal.
- Do not use kit past expiration date and do not intermix components from kits with different serial numbers.
- The SNAP device must be in a horizontal position on a flat surface while performing test.
- Do not use a SNAP device that has been activated prior to the addition of sample.
- WARNING:** Conjugate - H316/P332+P313/EUH208. Causes mild skin irritation. If skin irritation occurs: Get medical advice/attention. Contains Proclin. May produce an allergic reaction.

Storage

SNAP devices and test reagents must be stored at 2–8°C.

Kit components

Item	Reagents	Amount 15 Test	Amount 30 Test
1	bottle Anti-FeLV: HRPO Conjugate (Preserved with gentamicin and Proclin)	3.5 mL	7.0 mL
2	SNAP Device - Monoclonal antibodies to p27 and positive and negative controls	15	30
Reagents contained in each device:			
Wash Solution (preserved with Proclin)	0.4 mL	0.4 mL	
Substrate Solution	0.6 mL	0.6 mL	
Other components: transfer pipettes, sample tubes, and reagent rack			

Sample information

- Samples must be at room temperature (18–25°C) before beginning test procedure.
- Serum or plasma, either fresh, previously frozen or stored at 2–8°C, may be used in this test.
- Serum or plasma may be stored for up to 7 days at 2–8°C. For longer storage, sample should be frozen (20°C or colder).
- Previously frozen or older samples must be centrifuged before use.
- Hemolyzed or lipemic samples will not affect results.
- EDTA or heparin in plasma will not affect results.
- Whole blood may be used. Whole blood must be anticoagulated (eg. EDTA, heparin) and may be used either fresh or after refrigeration (2–8°C) for up to one week.

Test procedure

- Allow components to equilibrate at room temperature (18–25°C) for 30 minutes. **DO NOT HEAT.**

2. Using the pipet that is provided, transfer **3 drops of sample** (whole blood, serum, or plasma) into the sample tube.


3. Holding bottle vertical, add **4 drops of conjugate** to the sample tube.


4. Cap the sample tube and mix thoroughly by inverting tube **5-6 times**.

5. Place the device on a flat surface. Add contents of sample tube to sample well, being careful not to splash contents outside of sample well. Sample will flow across the result window, reaching the activation circle in approximately 30-60 seconds. Some sample may remain in sample well. **Watch carefully for sample or blue color in the activation circle.**

6. When fluid FIRST appears in the activation circle push activator firmly until it is flush with the device body.

Note: Some samples may not flow to the activation circle within 60 seconds, and, therefore, the circle may not turn color. In this case, press the activator after sample has flowed across result window. Keep the device horizontal to ensure accurate results.

7. Read test result at **10 minutes**.

NOTE: Positive control may develop sooner, but results are not complete until 10 minutes.

Interpreting test result

To determine test results, read the reaction spots in the result window. Any color development in the sample spots indicates the presence of FeLV antigen in the sample. If no color develops in the positive control spot, repeat the test.

Negative result

Only positive control spot develops color.

Positive result

If color in the sample spot is darker than negative control spot, result is positive.

Invalid results

*** Background**
If the sample is allowed to flow past the activation circle, background color may result. Some background color is normal. However, if colored background obscures test result, repeat the test.

*** No color development**
If positive control does not develop color, repeat the test.

Comparison	Sample size	SNAP test/Reference test	Sample type	Relative sensitivity and specificity 95% confidence limit (CL)	Kappa statistic		
PetChek® FeLV Ag	50	0	1	132 183	Serum/plasma/ whole blood	Spec., 100% (95% CL 91.3%–100%) Spec., 99.2% (95% CL 95.4%–100%)	0.99

Symbol descriptions

[LOT]	Batch code (lot)	[ECREP]	Authorized representative in the European Community
[SN]	Serial number	[IVD]	Consult instructions for use
[I]	Temperature limitation	[IVD]	In vitro diagnostic
[REF]	Catalog number	[REF]	Date of manufacture

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© 2018 IDEXX Laboratories, Inc. Alle rechten voorbehouden.