FIV Ab Single Test 10 x #(2)800005
Rev.: KI-40-00210-01

FIV Ab Single Test

Rapid Test

- (PT) Teste para FIV Ab
- (CS) FIV Ab Single Test
- (PL) Test FIV Ab
- (RO) Testul unic FIV Ab
- (SK) FIV Ab Single Test
- (DA) FIV Ab Single Test
- (NO) FIV Ab enkelt test
- (SV) Singeltest FIV Ab
- (ZH) FIV Ab 单项测试

FIV Ab Single Test

Rapid Test



- (ES) Test de único Ac de VIF
- (FR) FIV Ab Single Test
- (DE) FIV Ab Einzeltest
- (IT) Test singolo FIV Ab
- (NL) FIV Ab enkele test

For use in cats

Reorder # (2)800005


LOT 40002TTMMJJ

XXXX-MM

8 72017139957
Distributed by:
Covetrus BV
Beversestraat 23
5431 SL Cuijk (NL)
cbproducts@covetrus.com
Distrivel, S.L.
Registro entidad:
importadora HCMR-0097

FIV Ab Single Test 2 x #(2)800015
Rev.: KI-40-00202-01

FIV Ab Single Test

Rapid Test

- (ES) Test de único Ac de VIF
- (FR) FIV Ab Single Test
- (DE) FIV Ab Einzeltest
- (IT) Test singolo FIV Ab
- (NL) FIV Ab enkele test
- (PT) Teste para FIV Ab
- (CS) FIV Ab Single Test
- (PL) Test FIV Ab
- (RO) Testul unic FIV Ab
- (SK) FIV Ab Single Test
- (DA) FIV Ab Single Test
- (NO) FIV Ab enkelt test
- (SV) Singeltest FIV Ab
- (ZH) FIV Ab 单项测试

For use in cats

Reorder # (2)800015

Rev.: PO-40-002-01

covetrus

FIV Ab Single Test single teststrip

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(EN) FIV Ab Single Test Instructions for use

The Covetrus FIV Ab rapid test is used to test for feline immunodeficiency virus (FIV) antibodies in serum, plasma or whole blood.

Handling of samples

Do not refrigerate samples at room temperature for extended periods. Serum and plasma can be stored between 2-8°C for up to three days. For longer storage, the samples should be frozen at below -20°C. Venous whole blood should be stored between 2-8°C, if the test is due to be carried out within two days from sampling. Do not freeze whole blood samples. Bring the samples to room temperature before testing. Frozen samples must be fully thawed before testing and mixed well. Samples should not be refrozen and thawed.

General instructions for the use of blood samples

Separate the serum by centrifugation of the sample from the blood, to avoid hemolysis. • Heparinized and EDTA blood can be used to extract plasma. • Only use clear rather than haemolysed samples. • Testing should take place immediately after sampling.

Test procedure

Figure 1: Testing with serum and plasma

a. Draw blood sample from cat and extract the serum or plasma. Take the test pipette and draw sample material up into pipette. Add one drop of the sample material to the sample well. Wait a few seconds, to allow the strip to fully absorb the sample material.

b. Then add another drop of reagent to the test strip. If the liquid does not begin to run up on the strip after a few seconds, add another drop of reagent to the sample well.

Figure 2: Testing with whole blood

a. Draw blood sample from cat. Take the test pipette and draw sample material up into pipette. Add one drop of the sample material to the sample well. Wait up to one minute, to allow the strip to fully absorb the sample material.

b. Then add another drop of reagent to the test strip. If the liquid does not begin to run up on the strip after a few seconds, add another drop of reagent to the sample well.

Analysis (Figures 3 and 4):

After 10 minutes, one or two red lines will appear in the reaction field, and the test can be analysed.

Positive test result (Figure 3):

Both test line and the control line appear in the reaction field.

Negative test result (Figure 4):

Only the control line appears in the reaction field.

Invalid:

If no control line is shown, the test is invalid.

NB: All test results that are read after 10 minutes are invalid.

Disposal

All materials used must be disposed of carefully, because the material could be infectious. Test materials and used cassettes must be put in plastic bags and then disposed of.

Please note:

For veterinary use only. • Only for professional use. • Use the cassette as soon as possible after opening. • Please follow the procedure required for sample quantity. • Do not apply any sample material to the reaction field. • Do not touch the reaction field. • Use only the fluid supplied. • Do not use after the expiry date. • Do not use the test cassette if the packaging is damaged. • Observe analysis times indicated. Store and transport the test only at the temperatures indicated.

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