Leishmania Ab Single Test

Rapid Test

- (PT) Teste para Leishmaniose Ac
- (CS) Leishmania AB Single Test
- (PL) Test Leishmania AB
- (RO) Testul unic Leishmania AB
- (SK) Leishmania AB Single Test
- (DA) Leishmania AB Single Test
- (NO) Leishmania AB enkelt test
- (SV) Singeltest av antikroppar mot leishmaniasis
- (ZH) 利什曼原虫抗体单项测试

Leishmania Ab Single Test

Rapid Test

- (ES) Test único de Ac de Leishmania
- (FR) Leishmania AB Single Test
- (DE) Leishmania AB Einzeltest
- (IT) Test singolo Leishmania AB
- (NL) Leishmania AB enkele test

For use in dogs

Reorder # (2)800013

LOT

40010TTMMJJ


2022-03

2°C - 30°C

8 720171 399872

Distributed by:
Covetrus BV
Beversstraat 23
5431 SL Cuijk (NL)
cbproducts@covetrus.com
DistriNet, S.L.
Registro entidad
importadora HCMLR-0097

Rev.: PO-40-010-01

Leishmania Ab Single Test
single teststrip

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
EN Leishmania AB Single Test

Instructions for use

The Covetrus Leishmania rapid test is used to test Leishmania antibodies in serum, plasma or whole blood.

Handling the sample

Do not freeze samples at room temperature for extended periods. Serum and plasma can be stored between 2-8°C for up to three days. For longer storage, the samples should be frozen at below -20°C. Venous whole blood should be stored between -2-8°C if the test is due to be carried out within two days from sampling. Do not freeze whole blood samples. Bring the samples to room temperature before testing. Frozen samples must be fully thawed before testing and mixed well. Samples should not be refrozen and thawed.

General instructions for the use of blood samples: Separate the serum as rapidly as possible from the blood, to avoid haemolysis. • Heparin blood and EDTA blood can be used to extract plasma. • Only use clear rather than haemolysed samples. • Testing should take place immediately after sampling.

Test procedure

Figure 1: Testing with serum and plasma

a. Draw blood sample from dog and extract the serum or plasma. Take the test pipette and draw sample material up into pipette. Add one drop of the sample material to the sample well. Wait a few seconds, to allow the fluid to fully absorb the sample well.

b. Then add two drops of reagent to the sample well. If the liquid does not begin to run up on the strip after a few seconds, add another drop of reagent to the sample well.

Figure 2: Testing with whole blood

a. Draw blood sample from dog. Take the test pipette and draw sample material up into pipette. Add one drop of the sample material to the sample well. Wait up to one minute, to allow the fluid to fully absorb the sample well.

b. Then add two drops of reagent to the sample well. The fluid should begin to run up the test strip. If the liquid does not begin to run up on the strip after a few seconds, add another drop of reagent to the sample well.

Analysis (Figures 3 and 4):

After 10 minutes, one or two red lines will appear in the reaction field, and the test can be analysed.

Positive test result (Figure 3):

Both the control and the test line appear in the reaction field.

NB: Even a faint test line must be interpreted as a positive test result.

Negative test result (Figure 4):

Only the control line appears in the reaction field.

Invalid:

If no control line is shown, the test is invalid.

NB: All test results that are ready after 10 minutes are invalid.

Disposal:

All materials used must be disposed of carefully, because the material could be infectious.

Test materials and test cassettes must be put in plastic bags and then disposed of.

Please note:

For veterinary use only. • Only for professional use. • Use the cassette as soon as possible after opening. • Please follow the procedure for required sample quantity. • Do not apply any sample material to the reaction field. • Do not touch the reaction field. • Use only the fluid supplied. • Do not use after the expiry date. • Do not use the test cassette if the packaging is damaged. • Observe analysis times indicated. Store and transport the test only at the temperatures indicated.

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